| BioStation 2a | a | 1 |
|---------------|---|--------|
| | | 1 |
| | | 1 |
| | | 1 |
| | | 2 |
| 1. | | 2 |
| 2 | | 2 |
| 3 | | 2 |
| 0. 4 | | 2 |
| 5 | | ے ح |
| 5. 6 | | ך ג |
| 0. | | ך כ |
| | | ך כ |
| שו | | ך ע |
| | | 4 |
| | | 4 |
| | | 4 |

BioStation 2a



1



2

3

1/4

.

| ■ ア |
|---------------|
| 1 ABA |

1.

| , , , , | | | | | |
|---------|--|---|---|---|--|
| | | , | , | , | |
| | | | | | |

.

.

2.

, ·

٥

통신

3.

4.

6.

ОК

.

.

,

.

3/4

. ^

•

•

키 메뉴 ID
근태 모드
Y 작업 코드
경비

iD
iD
iD

- https://kb.supremainc.com/knowledge/

 \sim

| _ | |
|---|--|
| > | • |
| | |
| | |
| | ^ |
| | |
| | / . |
| | From: https://kb.supremainc.com/knowledge/ - |
| | Permanent link: https://kb.supremainc.com/knowledge/doku.php?id=ko:biostation_2a_menu |
| | Last update: 2023/11/01 17:19 |

.