Table of Contents

/

How to configure Daylight Saving Time	1
BioStar 2 Server	
BioStar 2 Client	3
Device	4

System Configuration, BioStar 2

How to configure Daylight Saving Time

You should configure Daylight Saving Time in three parts: BioStar 2 Server, BioStar 2 Client and Device. To utilize Daylight Saving Time, you must configure the same Daylight Saving Time for both BioStar 2 Client and device.

Applied version: BioStar 2.6.0 or above Supported device:

Device	Version
BioStation 2	V1.5.0 or above
BioStation A2	V1.4.0 or above
CoreStation	V1.1.0 or above
BioEntry P2	V1.1.0 or above
BioStation L2	V1.3.0 or above
BioEntry N2	V1.0.0 or above
BioEntry W2	V1.2.0 or above
FaceStation 2	V1.1.0 or above

* Entry device which firmware version is V2.x is not supported

New devices that support this feature:

- FaceStation F2 V1.0.0 or above
- X-Station 2 V1.0.0 or above

BioStar 2 Server

1. Go to Setting > DAYLIGHT SAVING TIME.

2/5

BioStar 2	🛞 Settings 🔒 Port	(i) About (?) Help			
	Settings				
USER DEVICE	ACCOUNT	PREFERENCE	CARD	CARD FORMAT	SERVER
	TRIGGER & ACTION	SCHEDULE	ALERT	HTTPS	CLOUD
	R= IMAGE LOG	USB AGENT	FACE GROUP MATCHING	EQ AUDIT TRAIL	DAYLIGHT SAVING TIME

2. Create new Daylight Saving Time by clicking "+Add" button.

 Setting 		Name	Start Time	End Time			+ Add
	1	Daylight Saving Time 1	April 1st Monday 00:00	June Last Saturday 00:00	,	¥	
	2	Daylight Saving Time 2	September 1st Monday 00:00	November Last Sunday 00:00	1	Ψ	
D	ST which	h is currently applied to the client cannot be edited or deleted.					
						_	

Add Daylight Sa	aving Time	×
 Name Start Date/Time End Date/Time 	Month Week Month Week Day of Week 00:00 Day of Week 00:00	
	Add Close	

- The maximum number of Daylight Saving Time you can create: Up to 200
- The maximum length of DST name: 48 letters

BioStar 2 Client

1. Go to **Setting > PREFERENCE**.

BioStar 2	2 💮 Settings 🕂 Port	t 🕧 About 🥐 Help			
DASH BOARD	Settings				
	Q, account	PREFERENCE	CARD	CARD FORMAT	SERVER
ELEVATOR O ZONE CONTROL	TRIGGER & ACTION	SCHEDULE	ALERT	HTTPS	CLOUD
	R= IMAGE LOG	USB AGENT	FACE GROUP MATCHING	AUDIT TRAIL	DAYLIGHT SAVING TIME

2. Select Daylight Saving Time you want to configure. You can choose maximum two periods for Daylight Saving Time. **Please note that you cannot configure two Daylight Saving Time which period are overlapped.**

← Preference							
Language / Time Zone							
• Language	English	•					
Time Zone	(UTC+9:00) Seoul, Tokyo, Osaka, Sapporo, Yakutsk	•					
Daylight Saving Time	Daylight Saving Time 1 (April 1st Monday 00:00 ~ June Last Saturday 00:00)	-					

- All users who login to BioStar 2 can configure Daylight Saving Time in Setting > Preference. In other word, a user without the privilege for monitoring can configure it as well.
- If Daylight Saving Time is configured in BioStar 2 Client (Setting > PREFERENCE), it cannot be modified or deleted. If you go to Setting > DAYLIGHT SAVING TIME, selected DST's icon is deactivated.

3. During the Daylight Saving Time, the time of event logs (Monitoring & Real-time log) will apply +1 hour considering DST and Timezone.

Device

If you want to configure DST in BioStar 2..

1. Go to the Device setting in BioStar 2.

2. Select Daylight Saving Time. You can select maximum two DST. Please note that you cannot configure two Daylight Saving Time which period are overlapped.

BioStar 2	🔅 🛞 Settings 📮 Port 🕧 Abou	t 🕐 Help		
DASH BOARD	← BioStation A2 54	531003 (192.168.14.207)		
A USER	Information			
R	• Name	BioStation A2 541531003 (192.168.14.207)	Group	Device Group B
DEVICE	Device ID	541531003	Device Type	BioStation A2
DOOR	Firmware Version	1.5.0 [2018/03/04 14:07:06] † Firmware Upgrade	Product Name	BSA2-OIPW
D	Kernel Version	1.3.0 [2018/03/04 14:08:17]	Hardware Version	1.0.0
ELEVATOR	Restore to default	All Without Network	Locked	Unlock
\odot	Time Zone	(UTC+9:00) Seoul, Tokyo, Osaka, Sapporo, Yakutsk 💌	✓ Time Synchronizati	on with Server
ZONE	 Daylight Saving Time 	Daylight Saving Time 1 (April 1st Monday 00:00 ~ June Last Saturday 00:00) 💌		
ACCESS CONTROL	System			
	Display Date	2018/04/02 III 16:40:38 x	Get Time	Set Time

- If Daylight Saving Time is configured in device menu, it can be modified or deleted unlike BioStar2 Client. If you go to Setting > DAYLIGHT SAVING TIME, selected DST's icon is activated.
- Master device (or Default device) can select maximum two Daylight Saving Time.
- Slave device cannot select Daylight Saving Time by itself but its time will be displayed based on Daylight Saving Time which is configured in Master Device.
- The Event log is managed based on UTC-0. However, the information on UTC, Time Zone and whether DST (Daylight Saving Time) is applied or not is managed together.

If you want to configure DST on the device directly..

1. Go to the Device menu on the device.

2. Go to **DEVICE > Daylight Saving Time**.

3. If you create new Daylight Saving Time, the DST will be automatically named as DST 1(Daylight Saving Time 1), DST 2(Daylight Saving Time 2) in turn. The name of DST is different depending on device model.

4. If you delete DST 1(Daylight Saving Time 1), DST 2(Daylight Saving Time 2) will be changed into DST 1(Daylight Saving Time 1).

If you reset the device without network setting, Daylight Saving Time will be remained without initialization.

* Visitor's start & expiry time depends on the server PC's Windows time zone. * Visitor's update follows the BioStar 2 time zone, DST applied. * Visitor's authentication time follows the device's time, which follows the BioStar 2 time zone even with DST applied.

From:

http://kb.supremainc.com/knowledge/ -

Permanent link: http://kb.supremainc.com/knowledge/doku.php?id=en:how_to_configure_daylight_saving_time&rev=1639978482

Last update: 2021/12/20 14:34