Table of Contents

How To Issue BioStar 2 QR	1
How to issue BioStar 2 QR	1
Step 1: Set Email Setting in BioStar 2 and enable [QR] option to send the QR code to user	
emails	1
Step 2: Add a X-Station 2 QR device to BioStar 2	3
Step 3: Go to USER menu and add QR	3
Step 4: Check your inbox. Download the attached QR in the email content, or print the QR	
	5

/

How To Issue BioStar 2 QR

Affect Products & Versions	
Device	X-Station 2 v1.0.1_20210329 or above
BioStar 2	BioStar v2.8.11 or above

From BioStar 2.8.11, QR/Barcode is available to be used as BioStar 2 User credentials.

X-Station 2 QR models(XS2-QDPB and XS2-QAPB) can read QR or Barcode and authenticate the user. The server can send an email with QR information to users and the users can download QR codes on their mobile or PC.

Users can have BioStar 2 QR code or even third-party QR code/barcode data and use them as their user credential.

- This document covers how to issue BioStar 2 QR.
- BioStar 2 QR is only available with X-Station 2 QR products.

For more details about X-Station 2, please visit [BioStar 2] New Product Introduction: X-Station 2.

How to issue BioStar 2 QR

You can generate QR codes in BioStar 2 system and use them as user credentials and below steps are the procedures for issuing the QR.

- Step 1: Set Email Setting in BioStar 2 and enable [QR] option to send the QR codes to user emails
- Step 2: Add X-Station 2 QR device to BioStar 2
- Step 3: Go to USER menu and add a new user, update BioStar 2 QR information with user email
- Step 4: Confirm the email, download the attached QR in the email, or print out the QR

Step 1: Set Email Setting in BioStar 2 and enable [QR] option to send the QR code to user emails

Direction: Settings>EMAIL SETTING>Email Contents Setting

• Email Contents Setting is for Visual Face Enrollment and QR Enrollment only.

Click [SMTP Setting].

		`
SMTP Setting	SMTP setting	
Test Mail Recipient Address	connect@suprema	
	• SMTP Setting • Test Mail Recipient Address	SMTP Setting SMTP setting Test Mail Recipient Address Connect@suprema Send Email

Refer to the SMTP Setting used by your email system and update the information correctly.

- Please check Security Type and update the port information in accordance with the security type.
- User name and Password must be matched.

ATP Option	
- Sender Information -	
SMTP Server Name	Suprema test
Description	This is only for test purpose
Server Address	Enter SMTP Server Address
Port(default:25)	Enter SMTP Server Port
• User Name	Kim@suprema.co.kr
 Password 	••••••
Security Type	TLS
Sender	supremaHR@suprema.co.kr



Enable QR and fill in the other contents.

OR

• QR	Use
Email Title	[Suprema] Your QR code has been generated
 Company Name 	Suprema Inc.
 Company Logo 	Upload Logo
- Contact	HR Team

Step 2: Add a X-Station 2 QR device to BioStar 2

XS2-QDPB and XS2-QAPB can scan QR or Barcode. Please check the product name and add the device to BioStar 2.8.11 or above.

BioSta	r 2 🛞 Settings 器 Port 🕜 About 🕐 Help
DASH BOARD	SEARCH DEVICE
8	
USER	E 🔜 All Devices
	BioLite N2 538747164 (192.168.13.141)
DEVICE	FaceStation 2 270 (192.168.13.14 4)
DOOR	FaceStation F2 547836017 (192.168.13.110)
	X-Station 2 547837141 (192.168.13.179)

Step 3: Go to USER menu and add QR

To generate the BioStar 2 QR code, go to USER. Add a new user and add an Email to the user.

2021/11/22 13:49		4/7		How To Issue BioStar 2 QR
← Add New Use	er			
Information				
	• Name	Helen	• Email	helenjang@supsrema.co.kr
8	DepartmentID	335	Title Telephone	
+ Add Photo	• Group	All Users	• Status	Active
	 Period Operator Level 	2001/01/01 00:00 ~	2030/12/31 23:59 III • Access Group	-
Click [QR/Barcode]	in Credential.			
Credential				
• 🗌 PIN				
- Auth Mode	Device Default			
Credential	+ Fingerprint 🍈	+ Face 🙂 + Visual Fac	* 🖸 + Card 🖬 +	Mobile + QR/Barcode

Select BioStar 2 QR and click [Enroll].

• You can select Input Type as random card ID or manual ID input.

roll Card		
QR/Barcode Type	BioStar 2 QR	•
Information		
Card ID	162346636437750645	Use User ID
• Input Type	Use random card ID	

Enroll	Cancel	
--------	--------	--

Click [Apply] to save the user data.

• Now the user receives the email from BioStar 2 and BioStar 2 synchronizes the user data to the connected devices.

⁻ http://kb.supremainc.com/knowledge/

Check the success message.



Step 4: Check your inbox. Download the attached QR in the email content, or print the QR

The user can check his or her email on mobile or PC. The attached QR is now available to be used.

[Suprema] You	r QR code	has been	generated
---------------	-----------	----------	-----------



SUP	emə

Here is your QR Code

Dear @suprema.co.kr,

You have been approved to access the site of Suprema Inc. using a QR code. Please download the QR code and scan on the device as guided below.



Scan QR Code



Please align the QR code to the guide beam at a 5~20cm distance from the sensor

Please inform us of any problems regarding the QR code. Contact: HR Team

For BioStar 2 QR, the device Operation mode is the same as Card Authentication mode.

If the user cannot find the QR in his or her email inbox, he or she can click [Resend] to re-send the email again.

 1:1 Security Level 	Security Level Device Default					
	Туре	Card Data Format	Sun	nmary		
	BioStar 2 QR	BioStar 2 QR	ID: 16236251	(Issue Count: 1)	Resend Block	

From: http://kb.supremainc.com/knowledge/ -Permanent link: http://kb.supremainc.com/knowledge/doku.php?id=en:how_to_issue_biostar_2_qr Last update: 2021/11/17 10:54

- http://kb.supremainc.com/knowledge/