


Table of Contents

BS2_StartMonitoringLogEx	1
Declaration	1
Parameter	1
Return Value	1
See Also	1

 **Fix Me!** This page is not fully translated, yet. Please help completing the translation.
(remove this paragraph once the translation is finished)

[Log Management API](#) > [BS2_StartMonitoringLogEx](#)

BS2_StartMonitoringLogEx

[+ V2.7.1] Starts Real-time log streaming including the temperature information.

Declaration

```
#include "BS_API.h"

int BS2_StartMonitoringLogEx(void* context, uint32_t deviceId,
OnLogReceivedEx ptrLogReceivedEx);
```

Parameter

- [In] *context* : Context
- [In] *deviceId* : Device ID
- [In] *ptrLogReceivedEx* : Callback function that is called when a new log is received

Return Value

If successfully done, `BS_SDK_SUCCESS` will be returned.
If there is an error, the corresponding error code will be returned.

See Also

[BS2_StopMonitoringLog](#)

From:
<https://kb.supremainc.com/kbtest/> - **BioStar Device SDK**

Permanent link:
https://kb.supremainc.com/kbtest/doku.php?id=en:bs2_startmonitoringlogex&rev=1599721518

Last update: **2020/09/10 16:05**